

Light lamps with magnifications

59271 Cold light lamp with magnification

- Light without shadow, it does not warm up.
- 3 Dioptres lent.
- Diameter: 125 mm.
- Arms balanced by springs.
- Adjustable head.
- Arms lenght: 1050 mm.
- Light power: 22 wat.
- Desk bracket.
- Lux at 30 cm: 2.000.



59284 Light lamp with magnification

- 3 Dioptres lent.
- Desk light lamp.
- Sliding protection cover.
- Light power: 9 wat.
- Lux at 30 cm: 900.



59275 Light lamp with magnification

- 3 Dioptres standard lent.
- Lent diameter: 100 mm.
- Arms balanced by springs.
- Adjustable head.
- Arms lenght: 900 mm.
- Light power: 40 wat.
- Desk bracket.
- Lux at 30 cm: 800.



59277 Manual light lamp with magnification

- Cold catode fluorescent light source.
- Long life (10.000 hours) and great brightness.
- 3 dioptres pure glass.
- Lux at 15 cm: 1.000.
- Battery (4x1,5V) operation or mains operation (converter included).



Germicidal UV lamps

- High efficiency filtered UV lamps provided in three different wavelengths: (254, 312 y 365 nm) or combined.

Applications:

254 nm: Fluorochemical, pesticides analysis, Mycology, proteins, steroids, mineralogy, archaeology, toxicology, mutation studies vitamins.

312 nm: Electrophoresis DNA/RNA.

365 nm: Fluorochemical, pesticides analysis, bacteria identification titration, Cultures, TLC, lipids, quality control, food pollution, archaeology, Mineralogy.

Model	Type	Tubes	nm	$\mu\text{W}/\text{CM}^2$	Price
58340	VL-4L	1 x 4 W	365	400	
58341	VL-4C	1 x 4 W	254	340	
58342	VL-4LC	1 x 4 W	365	350	
		1 x 4 W	254	265	
58360	VL-6L	1 x 6 W	365	700	
58361	VL-6C	1 x 6 W	254	710	
58362	VL-6LC	1 x 6 W	365	610	
		1 x 6 W	254	400	
58363	VL-6M	1 x 6 W	312	680	
58364	VL-6MC	1 x 6 W	312	580	
		1 x 6 W	254	400	
58370	VL-115L	1 x 15W	365	1100	
58371	VL-115C	1 x 15W	254	730	
58372	VL-115M	1 x 15W	312	1000	

More configurations 2x15W,... Consult



UV lamp stands

58386 Desk stand for VL-4 and VL-6 lamp.

58387 Desk stand for VL-115 lamp.

58388 Hand stand for VL-4 and VL-6 lamp.



TLC dark room

- Welcomes one or two hand held UV lamps (VL-6) in any of the three possible wavelengths.
- UV absorber shield to protect the user from the UV light.
- Lamps not included.

58384 Dark room for TLC

- To adapt up to 2 lamps VL-6 (Lamps not included)

