



**59110 Colony counter "SENSOR"**

- Electronic recording through any kind of marker.
- Pressure sensor system provides uniform sensitivity over the whole of the working surface.
- 4 digits LED display with zero reset button.
- Accept 10 cm to 15 cm petri dishes using a flexible adjuster.
- 1,5x lent, in an adjustable focusing rotatory arm 360°.
- Background plate can be changed to white or black for easy counting.
- Marker included.
- Dimensions: 30 x 33 x 10 cm.
- Weight: 5 kg.
- 230V-50 Hz.

*(Compleat equipment with auxiliary magnifying glass, flexible arm, marker, contrast plate)*

---

**Spare parts:**

59112 Spare lamp

59113 Marker

59114 1,5x lent

59115 Contrast plate

59116 Flexible arem

---

### Manual differential counters

- 5 Units and 8 units differential counters with totalizer.
- For simplifying differential blood cell counts.
- Representing the normal Schilling Group of Leucocytes: BASO, EOISN, SEG, LYMPH, MONO (70015) and MYELO, JUVEN, STAB (70017).
- Display from 0 up to 999 in each window for each key.
- A bell rings each time the totalizer records to 100 cells (100,200,...).
- Lateral CLEAR key (rotate handle).
- Weight: 950 gr (70015) and 1.250 gr (70017).



**70015** Differential counter 5 keys

**70017** Differential counter 8 keys

### 4 digits counter

- Manual counter with 4 digits, up to 999.
- Zero with a rotate handle.
- Weight: 35 gr.



**70010** 4 digits counter

### Electronic differential counter

- 8 units differential counter with totalizer.
- Representing the normal Schilling Group of Leucocytes: BASO, EOISN, SEG, LYMPH, MONO, MYELO, JUVEN, STAB.
- Every keystroke add 1 unit in the corresponding window and in the totalizer.
- An acoustic signal sounds when the totalizer records to 100 (100,200...)
- Corrective function with switch.
- Zero function.
- Lock function.
- Weight: 1.300 gr.
- 220V-50 Hz.



**70018** Electronic differential counter