



## Plastic tube racks

- Alphanumerical graduation.
- Dimensions: 235x110x75 mm.

Reference	Capacity N° tubes	Diameter máx.	Price Unit
100481	50	12 mm	
100482	50	16 mm	
100483	50	18 mm	
100484	18	20 mm	
100485	18	25 mm	
100486	18	30 mm	
100487	8	35 mm	
100488	8	40 mm	

## Plastic tube racks

- Racks made of polipropylene, easy – mounting, sterilized, can be used in water baths,...
- Assorted colours: White, blue, yellow, red.

Reference	Capacity N° tubes	Diameter máx.	Price Unit
100491	90	13 mm	
100492	60	16 mm	
100493	40	20 mm	
100496	24	25 mm	
100497	21	30 mm	



## Reversible tube racks - plastic

- Polipropylene tube racks, can be used for freeze.
- Alphanumeric graduation, it allows insert in one side 0,5/1,5/2ml microtubes and in the other side tubes and vials up to diameter 12mm. 48 or 80 positions models.

Reference	Description	Colour	Price
42145	Reversible tube rack for 48 tubes	Blue	
42146	Reversible tube rack for 48 tubes	Green	
42147	Reversible tube rack for 48 tubes	Red	
42148	Reversible tube rack for 48 tubes	Yellow	
42195	Reversible tube rack for 80 tubes	Blue	
42196	Reversible tube rack for 80 tubes	Green	
42197	Reversible tube rack for 80 tubes	Red	
42198	Reversible tube rack for 80 tubes	Yellow	




# Gradillas

## Tube racks

### Plastic tube racks with lids

- Polypropylene tube racks with lid.
- Autoclavables and can be used for freeze.
- Dimensions: 130x130 (Width x Depth).
- Assorted colours: natural, blue, green, pink, yellow,...

Reference	Description	Price
42280	Tube rack with lid for 81 tubes Up to 12,5mm Ø and 77mm height	
42285	Tube rack with lid for 64 tubes Up to 13,5mm Ø and 100mm height	
42286	Tube rack with lid for 49 tubes Up to 16,5mm Ø and 77mm height	
42287	Tube rack with lid for 49 tubes Up to 16,5mm Ø and 100mm height	
42288	Tube rack with lid for 25 tubes 15 ml "Falcon" type. Up to 17mm Ø and 126mm height	
42289	Tube rack with lid for 10 tubes 50 ml "Falcon" type. Up to 30mm Ø and 126mm height	



Racks to store tube racks in ultra-freezers

### Plastic microtube racks with lid

- Polypropylene tube rack 130x130x50 mm.
- Transparent lid with numbering.
- Capacity 81 microtubes up to 1,5 ml and cryovials up to 2 ml.
- Suitable for ultra-freezing -85°C.
- Bases in colours: Yellow, green, red, blue and natural.

Code	Dimensions mm	Capacity N° tubes	Volume tubes ml	Colour	Price
42290	130x130x50	81	1,5 – 2	Natural	
42291	130x130x50	81	1,5 – 2	Blue	
42292	130x130x50	81	1,5 – 2	Yellow	
42293	130x130x50	81	1,5 – 2	Green	
42294	130x130x50	81	1,5 – 2	Red	
42270	130x130x53	64	1,5 – 2	Natural	
42271	130x130x53	64	1,5 – 2	Blue	
42272	130x130x53	64	1,5 – 2	Yellow	
42273	130x130x53	64	1,5 – 2	Green	
42274	130x130x53	64	1,5 – 2	Red	



Racks to store tube racks in ultra-freezers

### Aluminium tube rack "Z" shape

- Tube racks made of anodized aluminium.
- Can be stored in freezers.
- Numerical identify.

Reference	Capacity N° tubes	Diameter Max.	Price
33012	50	12 mm	
33014	50	14,5 mm	
33017	50	17 mm	



## Aluminium racks for freeze

- Designed to store, carry and freeze samples.
- Made of anodized aluminium.


Reference	Model	Dimensions mm HxWxD	Capacity n° tubes	Diameter max. mm	Price
34040	LBL.4L.30	30 x 157 x 157	100	14	
34045	LBL.4L.45	45 x 157 x 157	100	14	
34033	LBL.3L.30	30 x 157 x 157	81	16	
34035	LBL.3L.45	45 x 157 x 157	81	16	

## Container –Racks for aluminium tube racks

- Storage and carriage system to assure an optimum use of space.
- Light structure made of aluminium.
- Wide range of models.

### Vertical containers (without tube racks)

Code	Model	Dimensions mm WxDxH	Capacity N° tube racks	Tube hight max	Price
30205	LBL.8G	170x162x575	8	55	
30210	LBL.6G	170x162x445	6	55	
30200	LBL.5G	170x162x540	5	80	
30220	LBL.4G	170x162x445	4	80	
30230	LBL.3G	170x162x445	3	120	
30240	LBL.4M	170x162x345	4	65	
30250	LBL.3M	170x162x345	3	90	
30260	LBL.3P	170x162x260	3	60	
30270	LBL.2P	170x162x260	2	95	



### Horizontal containers (without tube racks)

Code	Model	Dimensions mm WxDxH	Capacity N° tube racks	Tube hight max	Price
39200	LBL.12H	175x490x240	12	55	
39203	LBL.9G	175x490x240	9	75	
39205	LBL.6G	175x490x240	6	113	

## Containers-racks for freeze

- For storage in ultra-freezers of: microtubes 0,2 -0,5 -1,5,.. microplates, tube racks, falcon type tubes,....
- Light structure made of anodized aluminium.
- Suitable for plastic and cardboard tube racks.
- Frontal extraction.



## Microtubes containers

Code	Model	Dimensions mm W x D x H	Capacity n° tube racks	Tube rack máx. height	Price
39220	LBL.AG4	143x565x240	16	50	
39222	LBL. AP4	143x440x240	12	50	
39224	LBL.AG3	143x565x240	9	75	
39226	LBL.AP3	143x440x240	9	75	

## Microplate containers

- 39240** Microplate container LBL.APH  
Dimensions: 143x565x245.  
Capacity: 72 microplates.



## Containers for 15 and 50 ml tubes

Code	Model	Dimensions mm W x D x H	Capacity n° tubes	Price
39260	LBL.AFP	120x560x127	88 x 15 ml	
39262	LBL.AFG	120x560x127	42 x 50 ml	



## Vertical containers – racks for freeze

- Wide range of vertical containers for plastic, cardboard tube racks, microtiter plates, PCR plates, microtube racks....
- Made to order. (Consult).



## Microtube racks and racks for freeze

### Plastic microtube racks

- Autoclavables.
- Suitable for ultra-freezers.

Code	Dimensions mm	Capacity N° tubes	Volume tube ml	Colour	Price
42290	130x130x50	81	1,5 – 2	Natural	
42291	130x130x50	81	1,5 – 2	Blue	
42292	130x130x50	81	1,5 – 2	Yellow	
42293	130x130x50	81	1,5 – 2	Green	
42294	130x130x50	81	1,5 – 2	Red	
42250	130x130x32	144	0,2	Natural	
42251	130x130x32	144	0,2	Blue	
42252	130x130x32	144 </td <td>0,2</td> <td>Yellow</td> <td></td>	0,2	Yellow	
42253	130x130x32	144	0,2	Green	
42254	130x130x32	144	0,2	Pink	
42260	130x130x37	100	0,5	Natural	
42261	130x130x37	100	0,5	Blue	
42262	130x130x37	100	0,5	Yellow	
42263	130x130x37	100	0,5	Green	
42264	130x130x37	100	0,5	Pink	

## File cabinet – racks for tube racks

### Horizontal file cabinets


Code	Dimensions mm W x D x H	N° tube racks Capacity	Tubes ml	Price
30720	139 x 414 x 224	12	1,5/2	
30722	139 x 549 x 224	16	1,5/2	
30723	139 x 549 x 280	20	1,5/2	
30726	139 x 685 x 224	20	1,5/2	
30730	139 x 278 x 245	14	0,2	
30732	139 x 414 x 245	21	0,2	
30734	139 x 549 x 245	28	0,2	
30736	139 x 685 x 245	35	0,2	
30740	139 x 278 x 241	12	0,5	
30742	139 x 414 x 241	18	0,5	
30744	139 x 549 x 241	24	0,5	
30746	139 x 685 x 241	30	0,5	

### Vertical file cabinets

Code	Dimensions mm W x D x H	N° tube racks Capacity	Tubes ml	Price
30710	138 x 138 x 337	6	1,5/2	
30712	138 x 138 x 613	11	1,5/2	
30716	138 x 138 x 665	12	1,5/2	
30709	138 x 138 x 351	10	0,2	
30711	138 x 138 x 420	12	0,2	
30715	138 x 138 x 558	16	0,2	
30719	138 x 138 x 283	7	0,5	
30721	138 x 138 x 362	9	0,5	
30725	138 x 138 x 558	13	0,5	

### 30727 Vertical file cabinets

- With 11 adjustable height shelves.
- Dimensions: 140x140x610 mm.





### Microtube racks and racks for freeze

#### Cartboard microtube racks

- Suitable for Ultra-freezers.
- Different colors under request.


Code	Dimensions mm	Capacity N° tubes	Volume tubes ml	Colour	Price
42331	136x136x25	256	0,2	White	
42333	133x133x50	100	1,5 – 2	White	
42334	133x133x75	81	5	White	
42335	133x133x100	81	5	White	
42337	122x122x128	25	15	White	
42339	122x122x128	9	50	White	



### File cabinets for cartboard microtube racks

#### Horizontal file cabinets

Code	Dimensions mm WxDxH	N° microtube racks capacity	Tubes ml	Price
30760	139 x 414 x 224	12	1,5/2	
30762	139 x 549 x 224	16	1,5/2	
30763	139 x 549 x 280	20	1,5/2	
30765	139 x 549 x 224	12	5	
30770	142 x 425 x 213	21	0,2	
30772	142 x 565 x 243	32	0,2	
30774	142 x 565 x 273	36	0,2	
30780	127 x 380 x 265	6	15-50	
30782	127 x 505 x 265	8	15-50	



#### Vertical file cabinets

Code	Dimensions mm WxDxH	N° microtube racks capacity	Tubes ml	Price
30710	138 x 138 x 337	6	1,5/2	
30712	138 x 138 x 613	11	1,5/2	
30716	138 x 138 x 665	12	1,5/2	
30766	140 x 140 x 425	4	5	
30768	140 x 140 x 635	6	5	
30776	141 x 141 x 365	12	0,2	
30777	141 x 141 x 515	17	0,2	
30779	141 x 141 x 605	20	0,2	
30786	127 x 127 x 395	3	15-50	
30787	127 x 127 x 525	4	15-50	
30788	127 x 127 x 655	5	15-50	



#### 30727 Vertical file cabinet

- With 11 adjustable weight shelves.
- Dimensions: 140x140x610 mm.

30728 Additional shelf for file cabinet