



LABOLAN's laboratory desing department, is available to design your laboratory in accordance with your particular needs.

Our manufacturing, Installation and assembling department will culminate the design of your laboratory.

*INDELAB forniture basical characteristics, manufactured by LABOLAN:*

- Modular, easy installation, assembling and/or dis assembling.
- Metallic structure made of laminated stainless steel, protected by Epoxy paint covering, levelling elements lodges in its base and back side registrable.
- Worktop made of diverse materials with chemical resistance: granite, sandstone, resins, stratified, tempered glass, polypropilene, HPL.
- Modular sinks manufactured with the metallic structure described above.
- Modules under surface: suspended/movable with wheels or floor molding.
- Showcase cabinets and suspended cabinets.
- Safety cabinets for inflammables, acids and lyes according to norms in use.
- Laboratory gases aspiration and extraction cabinet.



*Ask us for our Modular Furniture Catalogue INDELAB*

### Standar laboratory bench

- Structure made of laminated stainless steel 60x30x1,5 mm, protected by epoxy paint covering.
- Back side registrable, levelling elements lodges in its base.
- Working surfaces:
  - LP: Laminated 0,8 mm over agglomerate 30 mm.
  - HPL: High pressure laminates compact 18 mm thickness.
  - TPL: Trespa Top Lab Plus 20 mm thickness (High resistance to high concentration acids pigments,...)
  - (Other working surfaces: Stainless steel, granite, tempered glass, Vitrified monolithic sandstone,.....consult).

Reference	Dimensions cm WxDxH	Worktop	Price
97120750	120 x 75 x 90	LP	
98120750	120 x 75 x 90	HPL	
99120750	120 x 75 x 90	TLP	
97150750	150 x 75 x 90	LP	
98150750	150 x 75 x 90	HPL	
99150750	150 x 75 x 90	TLP	
97180750	180 x 75 x 90	LP	
98180750	180 x 75 x 90	HPL	
99180750	180 x 75 x 90	TLP	

### Housing for services

- Made of rolled steel protected by epoxy paint covering.
- Double adjustable height shelf.
- Gallery of services.
- Electrical conduits with bases (optional).

- 96120015 Mural 120x15 cm.  
 96150015 Mural 150x15 cm.  
 96180015 Mural 180x15 cm.  
 96120030 Central 120x30 cm.  
 96150030 Central 150x30 cm.  
 96180030 Central 180x30 cm.

### Under bench furniture

- Modules under bench made of melamina 19 mm.
- Edges made of PVC 3 mm thickness.
- Drawers made of high resistance polipropylene.  
 Versions:
  - "C" Modules under surface suspended and movables.
  - "R" Modules under surface with wheels.
  - "Z" Modules under surfaces with floor molding.

Code	Description	Dimensions	Price
9460451 -	Modules under bench door/pane adjustable	60x45x65 cm	
9460452 -	Modules under bench drawer/door/pane adjustable	60x45x65 cm	
9460453 -	Modules under bench 4 drawers	60x45x65 cm	

(Replace " - " with the version of the module you want to order: "C"- "R" o "Z")

(Other dimensions or versions, consult)

### Washing/sink units

- Made of rolled steel protected by epoxy paint covering.
- Versions for sink and wiper:
  - Stainless steel, polipropylene acidproof and vitrified monolithic sandstone.



## Laboratory fittings

- According to DIN 12920, 12898, 12919 and 3537 norms.

**73111** Water tap throat, gyrotary, with 1 outlet

**73113** Water tap throat, gyrotary, with 2 outlets

**73119** Water tap throat, gyrotary, mixer cold/hot

**77094** Water jet Vacuum pump

- High chemical resistance.
- Tap water inlet 3/4" included reducer 1/2", connection olive hose.
- Vacuum connection unscrewed olive.
- Length: 210 mm.

**73220** Bench mounted Gas tap 1 outlet

**73222** Bench mounted Gas tap 2 outlets 180°

**73224** Bench mounted Gas tap 2 outlets 90°

**73320** Bench mounted vacuum tap 1 outlet

**73322** Bench mounted vacuum tap 2 outlets

**73901** Electrical turret with 2 outputs

**73902** Electrical turret with 4 outputs



## Washing units

- Washing units to be adapted in worktops or in independent modules.
- Available washing units made of different materials monolithic sandstone, stainless steel,...

**73301** Ceramic sink 15x15 cm

**73302** Ceramic sink 30x15 cm

**73400** Polypropylene sink 15x15 cm

**73415** Polypropylene sink 30x15 cm

**73420** Polypropylene washing unit 60x75 cm

- Sink dimensions: 50x40x32 cm.

**73425** Polypropylene washing unit 120x75 cm

- Sink dimensions 50x40x32 cm and wiper.



## Emergency Eye Wash

**72212** Emergency eye wash for laboratory

- Two high flow water sprays aerated at low pressure.
- The valve stays open until it is manually closed.

**74220** Emergency Shower/eye wash

- Twin heads-ABS shower head.



- 73600 Stool with 5 legs and wheels or bumper**
- Footrest ring and adjustable height.
  - Without back and gyrotory with gas elevation.
  - Upholster seat.

- 73601 Stool with 5 legs and wheels or bumper**
- Footrest ring and adjustable height.
  - Without back and gyrotory with gas elevation.
  - Black skay seat.

- 73610 Stool with 5 legs and wheels or bumper**
- Footrest ring and adjustable height.
  - With back and gyrotory with gas elevation.
  - Upholster seat.

- 73611 Stool with 5 legs and wheels or bumper**
- Footrest ring and adjustable height.
  - With back and gyrotory with gas elevation.
  - Black skay seat.

- 73700 Chair for low desk with wheels or bumper**
- With back and gyrotory with gas elevation.
  - Black skay seat.

- 73641 Stool with 5 legs and wheels or bumper**
- Footrest ring.
  - Without back and gyrotory with gas elevation.
  - Black skay seat or upholster seat.

- 73651 Stool with 5 legs and wheels or bumper**
- Footrest ring.
  - With back and gyrotory with gas elevation.
  - Black skay seat or upholster seat.

- 73501 Stool with 4 legs and footrest ring**
- Gyrotory, adjustable height.
  - Black skay seat.

- 73502 Stool with 4 legs and footrest ring**
- Gyrotory, adjustable height.
  - Black polypropylene seat.

- 73511 Stool with 5 legs with back, gyrotory**
- With back, gyrotory, adjustable height.
  - Black skay seat.

